

# Taney County Ambulance District

## Public Access Defibrillation Program

### AED Site Information Form

(Please complete one form for each unit)

---

**Company Name:** \_\_\_\_\_

**Address:** \_\_\_\_\_ **City:** \_\_\_\_\_ **State** \_\_\_\_\_

**Business Phone Number:** \_\_\_\_\_ **Fax #:** \_\_\_\_\_

**Type of Business:** \_\_\_\_\_

**AED Location:** \_\_\_\_\_ **In a Cabinet:** Yes No

**AED Contact Person/Program Coordinator:** \_\_\_\_\_

**Contact Person Phone Number:** \_\_\_\_\_ **Ext:** \_\_\_\_\_

**Contact Person E-mail:** \_\_\_\_\_

---

**Brand of AED:** \_\_\_\_\_ **Model #:** \_\_\_\_\_

**Date AED Installed:** \_\_\_\_\_ **Serial #:** \_\_\_\_\_

**Number of Trained Employees:** \_\_\_\_\_

**Circle Type(s) of Electrodes Located with AED:** Adult Child Ped

Please return this Form to: Taney County Ambulance District

106 Industrial Park Rd

Hollister, MO. 65672

417-339-2898 FAX – 417-336-3727

## **Taney County Ambulance District**

### **Ordinance 97-01 Ensuring Public Health and Welfare**

#### **Chapter 400 – AED (Automated External Defibrillator)**

##### 97-01.401 – Health and Safety Standards

**PROPOSE:** This rule provides the health and safety requirements and standards as related to AED's.

1. All AED's will be registered with Taney County Ambulance District.
2. All AED's should be maintained and serviced according to state and local laws.
3. All public and private entities must follow TCAD's or approved AED protocols for usage.
4. All vendors will be required to notify TCAD upon the sale of any AED's in Taney County.
5. Upon the use of any AED within Taney County, a report must be filed with tCAD within 24 hours of its use.
6. TCAD reserves the right to inspect any AED's within Taney County for the following:
  - a. Proper training of personnel.
  - b. Proper maintenance of unit.
  - c. Proper storage of the unit.
  - d. Up to date electrode patches.
  - e. Medical protocols kept with unit.
7. The Quality Insurance Committee will review all AED usage, during their regular scheduled meetings.

# Taney County Ambulance District

## Policy and Procedure Manual

<u>Chapter 300</u>	Patient Care Protocols
<u>Section 03</u>	Cardiac Emergencies
<u>Title</u>	<b>Automatic External Defibrillation</b>
<u>Number</u>	303-8

Distribution – All Operations Personnel

Purpose – To outline the AED treatment protocol for patients suffering from cardiac arrest.

### Protocol

1. Immediately upon arrival at the scene, assess the patient and verify cardiac arrest by the absence of consciousness, respiration and pulse.
2. Initiate CPR, ventilating with 100% Oxygen.
3. Prepare the patient while continuing with CPR:
  - a. Expose the patient's chest by removing clothing.
  - b. Dry the chest and remove nitroglycerin patches if necessary
  - c. If chest hair is excessive, quickly shave the area to permit proper adhesion of the electrodes.
4. Prepare the electrodes:
  - a. Remove the electrodes from the packaging.
  - b. Attach the defibrillator cables to the electrodes.
5. Attach electrodes to the patient; turn the Automated External Defibrillator (AED) power on. Follow prompts from the AED.
6. Follow the prompts of the AED making sure that no one is touching the patient or equipment during prompted defibrillations.

7. Follow current AHA guidelines for CPR. **REMEMBER, GOOD CPR IS CRUTIAL.**
8. If patient has return of spontaneous signs of circulation, refer to post resuscitation protocol (303-7).
9. Begin transport as soon as possible if available or notify the responding ALS unit and receiving E.D. of patient status/condition as soon as possible.

Medical Director Signature: Dr. L.E. Mire

Date: 10/28/08